

Product datasheet

Goat Anti-Alligator IgG H&L (HRP) ab112780

Overview

Product name	Goat Anti-Alligator IgG H&L (HRP)
Description	Goat polyclonal Secondary Antibody to Alligator IgG - H&L (HRP)
Host species	Goat
Target species	Alligator
Specificity	ab112780 reacts specifically with Alligator IgG and with light chains common to other Alligator immunoglobulins. No antibody was detected against nonimmunoglobulin serum proteins.
Tested applications	Suitable for: WB, ELISA, IHC-P, ICC/IF
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
Purity	Immunogen affinity purified
Purification notes	ab112780 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to horseradish peroxidase (HRP).
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab112780 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000 - 1/30000. Colorimetric: 1/1000-1/10000. Chemiluminescence: 1/1000-1/30000.

Application	Abreviews	Notes
ELISA		1/15000 - 1/50000. Coating.
IHC-P		1/200 - 1/500.
ICC/IF		1/200 - 1/500.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors