

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG - (Fab')₂ (HRP), pre-adsorbed ab98762

Overview

Product name	Goat F(ab') ₂ Anti-Mouse IgG - (Fab') ₂ (HRP), pre-adsorbed	
Host species	Goat	
Target species	Mouse	
Tested applications	Suitable for: WB, ELISA, IHC-P, ICC	
Minimal cross-reactivity	Human, Rat	more details
Conjugation	HRP	

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to Horse Radish Peroxidase.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab98762 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/5000 - 1/25000. by chemiluminescent detection. Use at 1/2000 - 1/20000 by colorimetric detection.

Application	Abreviews	Notes
ELISA		1/10000 - 1/100000. as primary
IHC-P		1/200 - 1/2000.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors