

Product datasheet

Goat Anti-Horse IgGT H&L (FITC) ab112874

Overview

Product name	Goat Anti-Horse IgGT H&L (FITC)
Description	Goat polyclonal Secondary Antibody to Horse IgGT (FITC)
Host species	Goat
Target species	Horse
Specificity	By immunoelectrophoresis and ELISA ab112874 reacts specifically with Horse IgG(T). It may cross react with IgG from other species.
Tested applications	Suitable for: ICC/IF, IHC-P
Immunogen	Horse IgG(T)
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 99% PBS
Purity	Immunogen affinity purified
Purification notes	Antiserum was solid phase adsorbed to ensure subclass specificity. ab112874 was isolated by affinity chromatography using antigen coupled to agarose beads.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab112874 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/200 - 1/2000.

Application	Abreviews	Notes
IHC-P		1/200 - 1/2000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors