

Product datasheet

Goat Anti-Human IgG Fc (HRP) preadsorbed ab205630

Overview

Product name	Goat Anti-Human IgG Fc (HRP) preadsorbed
Host species	Goat
Target species	Human
Tested applications	Suitable for: ELISA
Immunogen	Full length protein corresponding to Human IgG. normal human IgG and human IgG paraproteins
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot. Store at -20°C. Stable for 12 months at -20°C. Store in the Dark.
Storage buffer	pH: 7.40 Constituents: 50% Glycerol (glycerin, glycerine), 50% PBS
Purity	IgG fraction
Purification notes	Source: Pooled antisera from goats hyperimmunized with normal human IgG and human IgG paraproteins;
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab205630** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/4000 - 1/8000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors