

## Product datasheet

# Goat Anti-Llama IgG H&L (FITC) ab112785

### Overview

<b>Product name</b>	Goat Anti-Llama IgG H&L (FITC)
<b>Description</b>	Goat polyclonal Secondary Antibody to Llama IgG - H&L (FITC)
<b>Host species</b>	Goat
<b>Target species</b>	Llama
<b>Specificity</b>	By immunoelectrophoresis and ELISA ab112785 reacts specifically with llama IgG and with light chains common to other llama immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. ab112785 may cross react with IgG from other species.
<b>Tested applications</b>	<b>Suitable for:</b> ICC/IF, IHC-P
<b>Conjugation</b>	FITC. Ex: 493nm, Em: 528nm

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 98% PBS, 0.2% BSA
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	The antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab112785 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/200 - 1/2000.

Application	Abreviews	Notes
IHC-P		1/200 - 1/2000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors