

Product datasheet

Goat Anti-Mouse IgG Fc (18nm Gold) preadsorbed ab105287

Overview

Product name	Goat Anti-Mouse IgG Fc (18nm Gold) preadsorbed
Host species	Goat
Target species	Mouse
Specificity	Recognizes the heavy chains on mouse IgG. Exhibits minimal crossreactivity to human, bovine and horse serum proteins.
Tested applications	Suitable for: ELISA, Electron Microscopy, WB
Minimal cross-reactivity	Cow, Horse, Human more details
Conjugation	Gold 18nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 0.87% Sodium chloride, 0.164% Sodium phosphate, 0.381% Sodium borate, 1.5% BSA
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab105287** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.

Application	Abreviews	Notes
Electron Microscopy		1/10 - 1/20.
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors