

Product datasheet

Goat Anti-Rat IgG H&L (HRP polymer) ab214882

★★★★★ 3 Abreviews 3 References 1 Image

Overview

Product name	Goat Anti-Rat IgG H&L (HRP polymer)
Host species	Goat
Target species	Rat
Specificity	Based on IEP, ab214882 reacts with heavy (gamma) chains on rat IgG, light chains on all rat immunoglobulins. No cross-reactivity is observed to non-immunoglobulin rat serum proteins, IgG from human or mouse.
Tested applications	Suitable for: Immunoelectrophoresis, IHC-P
Immunogen	Full length protein corresponding to Rat IgG. (Purified Rat IgG, whole molecule).
Conjugation	HRP polymer

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.1% Proclin 300 Constituents: 1% BSA, 98% PBS
Purity	Affinity purified
Purification notes	ab214882 was affinity purified using solid phase Rat IgG. Purity >95% by SDS-PAGE.
Clonality	Polyclonal
Isotype	IgG


Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab214882 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes

Application	Abreviews	Notes
Immunoelectrophoresis		Use at an assay dependent concentration. This product is provided in a pre-diluted, ready-to-use format. Incubate 150 µl per slide for 1 hour.
IHC-P	★★★★★ (2)	Use at an assay dependent concentration.

Images



C57BL/6 mouse spleen stained with 1 µg/ml of rat anti-mouse CD45R/B220 detected using ab214882 Goat Anti-Rat IgG H&L (HRP polymer).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Rat IgG H&L (HRP polymer) (ab214882)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors