

## Product datasheet

# Goat Anti-Rhesus monkey IgG H&L (Rhodamine) ab125851

### Overview

<b>Product name</b>	Goat Anti-Rhesus monkey IgG H&L (Rhodamine)
<b>Host species</b>	Goat
<b>Target species</b>	Rhesus monkey
<b>Tested applications</b>	<b>Suitable for:</b> ICC/IF, ELISA, WB, Dot blot, IHC-P
<b>Immunogen</b>	Full length native protein (purified) corresponding to IgG. Pooled antisera from goats hyperimmunized with rhesus IgG
<b>Conjugation</b>	Rhodamine. Ex: 550nm, Em: 570nm

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: 99% PBS
<b>Purity</b>	Affinity purified
<b>Purification notes</b>	Affinity chromatography on rhesus IgG covalently linked to agarose.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab125851 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use 1µg for 10 <sup>6</sup> cells.
ELISA		Use at an assay dependent concentration.

Application	Abreviews	Notes
<b>WB</b>		Use at an assay dependent concentration.
<b>Dot blot</b>		Use at an assay dependent concentration.
<b>IHC-P</b>		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors