

Product datasheet

GOLGA3 peptide (1486-1498) (Carboxyterminal end)
 ab45813

1 Image

Overview

Product name GOLGA3 peptide (1486-1498) (Carboxyterminal end)

Description

Nature Synthetic

Amino Acid Sequence

Sequence HSQSRASKEGPGE

Amino acids 1486 to 1498

Specifications

Our [Abpromise guarantee](#) covers the use of **ab45813** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Blocking - Blocking peptide for Anti-GOLGA3 antibody ([ab40837](#))

Form Liquid

Preparation and Storage


Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

General Info

Relevance GOLGA3 is a member of the golgin family of proteins which are localized to the Golgi. It has been postulated to play a role in nuclear transport and Golgi apparatus localization.

Cellular localization Cytoplasmic and Golgi Apparatus

Images



Western blot - GOLGA3 peptide (1486-1498)
(Carboxyterminal end) (ab45813)

All lanes : GOLGA3 peptide (1486-1498) (Carboxyterminal end)
(ab45813) at 0.1 µg/ml

Lane 1 : HeLa lysate

Lane 2 : HeLa lysate with Immunising peptide

Lysates/proteins at 35 µg per lane.

Developed using the ECL technique.

Observed band size: 170 kDa

[why is the actual band size different from the predicted?](#)

Additional bands at: 25 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors