

Product datasheet

Hamster (Armenian) Complement Preserved Serum ab155107

Overview

| | |
|----------------------|--|
| Product name | Hamster (Armenian) Complement Preserved Serum |
| Description | Hamster (Armenian) Complement Preserved Serum |
| Host species | Armenian hamster |
| General notes | Hamster Complement serum is collected from spontaneously clotted whole blood and centrifuged to separate the serum from the clot. Hamster Complement is ideal for hemolytic assay where hemolysis of IgG antibodies are required. Complement is preservative free. Raised in: Armenian Hamster Purity: Whole antiserum |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| Purity | Whole antiserum |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors