

## Product datasheet

# Anti-hCG Holo c2 antibody [INN-hCG-53] ab11392

### 1 References

#### Overview

---

<b>Product name</b>	Anti-hCG Holo c2 antibody [INN-hCG-53]
<b>Description</b>	Mouse monoclonal [INN-hCG-53] to hCG Holo c2
<b>Host species</b>	Mouse
<b>Specificity</b>	Ab11392 recognises epitope c2 on hCG. It does not recognise free hCG subunits.
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	This information is proprietary to Abcam and/or its suppliers.
<b>General notes</b>	<p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: PBS</p>
<b>Purity</b>	Protein G purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	INN-hCG-53
<b>Isotype</b>	IgG

#### Target

---

<b>Relevance</b>	The placenta produces hCG, whose presence, along with other hormones, helps maintain the early stages of pregnancy. Because hCG is produced only by placental tissue and the hormone can be found in the blood or urine of a pregnant woman, it has become a convenient chemical test of pregnancy.
<b>Cellular localization</b>	Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors