

Product datasheet

Hexanoyl-Lysine BSA Conjugate ab245918

[1 Image](#)

Description

Product name	Hexanoyl-Lysine BSA Conjugate
Purity	> 98 % SDS-PAGE.
Expression system	Synthetic
Protein length	Full length protein
Animal free	No
Nature	Native
Species	Cow
Additional sequence information	Bovine Serum Albumin (BSA) modified with Hexanoyl-Lysine adduct (HEL).
Description	Native Cow Bovine Serum Albumin protein

Specifications

Our **[Abpromise guarantee](#)** covers the use of **ab245918** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Applications	Western blot ELISA SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle. Preservative: 0.09% Sodium azide Constituent: PBS Filter sterilized (0.2 µm)
------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Images



Western blot - Hexanoyl-Lysine BSA Conjugate (ab245918)

All lanes : Anti-Hexanoyl-Lysine adduct Antibody, Clone 5D9, 2 hrs at RT at 1/1000 dilution

Lane 1 : Hexanoyl-Lysine BSA Conjugate (ab245918) at 1 μ g

Lane 2 : Hexanoyl-Lysine BSA Conjugate (ab245918) at 0.5 μ g

Secondary

All lanes : Goat Anti-Mouse IgG: HRP, 60 min at RT at 1/2000 dilution

Block: 5% Skim Milk in 1X TBST.

Color Development: Luminol for 1 min at RT.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors