

## Product datasheet

# Anti-Histone H4 (acetyl K5 + K8 + K12) antibody ab233193

[5 Images](#)

### Overview

---

<b>Product name</b>	Anti-Histone H4 (acetyl K5 + K8 + K12) antibody
<b>Description</b>	Rabbit polyclonal to Histone H4 (acetyl K5 + K8 + K12)
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> ChIP-sequencing, WB, ICC/IF, Dot blot, ChIP
<b>Species reactivity</b>	<b>Reacts with:</b> Mouse, Human, Recombinant fragment
<b>Immunogen</b>	Synthetic peptide corresponding to Human Histone H4 (acetyl K5 + K8 + K12) conjugated to keyhole limpet haemocyanin. Database link: <a href="#">P62805</a>
<b>Positive control</b>	ChIP: Chromatin from K562 cells. ChIPseq: Chromatin from K562 cells. WB: HeLa whole cell and histone extracts. ICC/IF: HeLa cells.
<b>General notes</b>	<p>Reproducibility is key to advancing scientific discovery and accelerating scientists' next breakthrough.</p> <p>Abcam is leading the way with our range of recombinant antibodies, knockout-validated antibodies and knockout cell lines, all of which support improved reproducibility.</p> <p>We are also planning to innovate the way in which we present recommended applications and species on our product datasheets, so that only applications &amp; species that have been tested in our own labs, our suppliers or by selected trusted collaborators are covered by our Abpromise™ guarantee.</p> <p>In preparation for this, we have started to update the applications &amp; species that this product is Abpromise guaranteed for.</p> <p>We are also updating the applications &amp; species that this product has been “predicted to work with,” however this information is not covered by our Abpromise guarantee.</p> <p>Applications &amp; species from publications and Abreviews that have not been tested in our own labs or in those of our suppliers are not covered by the Abpromise guarantee.</p> <p>Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, as well as customer reviews and Q&amp;As.</p>

## Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservatives: 0.05% Sodium azide, 0.05% Proclin 300 Constituent: PBS
<b>Purity</b>	Affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

## Applications

Our [Abpromise guarantee](#) covers the use of **ab233193** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ChIP-sequencing		Use at an assay dependent concentration. Use 0.5 µg.
WB		1/1000.
ICC/IF		1/500.
Dot blot		1/20000.
ChIP		Use at an assay dependent concentration. Use 0.5 - 1 µg per IP.

## Target

<b>Function</b>	Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.
<b>Sequence similarities</b>	Belongs to the histone H4 family.
<b>Post-translational modifications</b>	Acetylation at Lys-6 (H4K5ac), Lys-9 (H4K8ac), Lys-13 (H4K12ac) and Lys-17 (H4K16ac) occurs in coding regions of the genome but not in heterochromatin. Citrullination at Arg-4 (H4R3ci) by PADI4 impairs methylation. Monomethylation and asymmetric dimethylation at Arg-4 (H4R3me1 and H4R3me2a, respectively) by PRMT1 favors acetylation at Lys-9 (H4K8ac) and Lys-13 (H4K12ac). Demethylation is performed by JMJD6. Symmetric dimethylation on Arg-4 (H4R3me2s) by the PRDM1/PRMT5 complex may play a crucial role in the germ-cell lineage. Monomethylated, dimethylated or trimethylated at Lys-21 (H4K20me1, H4K20me2, H4K20me3). Monomethylation is performed by SET8. Trimethylation is performed by SUV420H1 and SUV420H2 and induces gene silencing.


Ubiquitinated by the CUL4-DDB-RBX1 complex in response to ultraviolet irradiation. This may weaken the interaction between histones and DNA and facilitate DNA accessibility to repair proteins. Monoubiquitinated at Lys-92 of histone H4 (H4K91ub1) in response to DNA damage. The exact role of H4K91ub1 in DNA damage response is still unclear but it may function as a licensing signal for additional histone H4 post-translational modifications such as H4 Lys-21 methylation (H4K20me).

Sumoylated, which is associated with transcriptional repression.

## Cellular localization

Nucleus. Chromosome.

## Images



**All lanes :** Anti-Histone H4 (acetyl K5 + K8 + K12) antibody (ab233193) at 1/1000 dilution

**Lane 1 :** HeLa (human epithelial cell line from cervix adenocarcinoma) whole cell extract at 25 µg

**Lane 2 :** HeLa histone extract at 15 µg


**Lane 3 :** Recombinant histone H2A at 1 µg

**Lane 4 :** Recombinant histone H2B at 1 µg


**Lane 5 :** Recombinant histone H3 at 1 µg

**Lane 6 :** Recombinant histone H4 at 1 µg

**Dilution buffer:** TBS-Tween containing 5% skimmed milk.




ChIP assays were performed using human K562 (Human chronic myelogenous leukemia cell line from bone marrow) cells, ab233193 and optimized PCR primer sets for qPCR. ChIP was performed on sheared chromatin from 100,000 cells. A titration of the antibody consisting of 0.2, 0.5, 1 and 2 µg per ChIP experiment was analysed. IgG (1 µg/IP) was used as negative IP control. QPCR was performed with primers for promoter of the active gene EIF4A2 and for a region 1 kb upstream of the GAPDH gene, used as positive controls, and for the inactive MYOD1 gene and the Sat2 satellite repeat region used as negative controls. Image shows the recovery, expressed as a % of input (the relative amount of immunoprecipitated DNA compared to input DNA after qPCR analysis).



ChIP-sequencing - Anti-Histone H4 (acetyl K5 + K8 + K12) antibody (ab233193)

ChIP was performed with 0.5  $\mu$ g ab233193 on sheared chromatin from 100,000 K562 cells. The IP'd DNA was subsequently analysed on an Illumina Genome Analyzer. Library preparation, cluster generation and sequencing were performed according to the manufacturer's instructions. The 36 bp tags were aligned to the human genome using the ELAND algorithm. The figure shows the signal distribution along the complete length of chromosome 2 (figure A) and a zoomin to a 600 kb region (figure B). Figure C and D show the enrichment in two genomic regions on chromosome 3 and 12, respectively, containing EIF4A2 and GAPDH positive controls.




Immunocytochemistry/ Immunofluorescence - Anti-Histone H4 (acetyl K5 + K8 + K12) antibody (ab233193)

HeLa (Human epithelial cell line from cervix adenocarcinoma) cells were stained with ab233193 and with DAPI.

Cells were fixed with 4% formaldehyde for 10 minutes and blocked with PBS/TX-100 containing 5% normal goat serum and 1% BSA.

**Figure A:** Cells were immunofluorescently labeled with ab233193 (left) diluted 1/500 in blocking solution followed by an anti-rabbit antibody conjugated to Alexa<sup>®</sup>488. The middle panel shows staining of the nuclei with DAPI. A merge of the two stainings is shown on the right.

**Figure B, C, D and E:** Staining of the cells with ab233193 after incubation of the antibody with 10 ng/ $\mu$ l of the following blocking peptides: H4K5,8,12 unmodified (B), H4K5,8,12ac (C), H2A.ZK5,7,11ac (D) and H4K5,8,12,16ac (E).



To test the cross reactivity of ab233193 a Dot Blot analysis was performed with peptides containing other histone modifications and the unmodified H4. One hundred to 0.2 pmol of the respective peptides were spotted on a membrane. ab233193 was used at a dilution of 1:20,000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors