

Product datasheet

Anti-HLA Class II DRB1 antibody (FITC) ab193620

Overview

Product name	Anti-HLA Class II DRB1 antibody (FITC)
Description	Rabbit polyclonal to HLA Class II DRB1 (FITC)
Host species	Rabbit
Conjugation	FITC. Ex: 493nm, Em: 528nm
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment corresponding to Human HLA Class II DRB1 aa 31-266. Sequence: DTRPRFLEYSTSECHFFNGTERVRYLDRYFHNQEENV RFDSVGEFRAVT ELGRPDAEYWNSQKDLLEQKRGRVDNYCRHNYGVVE SFTVQRRVHPKVTV YPSKTQPLQHLLVCSVSGFYPGSIEVRWFRNGQEE KTVVSTGLIHNG DWTFTLVMLETVPRSGEVYTCQVEHPSVTSPLTVE WRARSESAQSKMLS GVGGFVLGLLFLGAGLFIYFRNQKGHSGLQPRGFLS

Database link: [P01912](#)

 [Run BLAST with](#)

 [Run BLAST with](#)

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.4 Preservative: 0.03% Proclin Constituents: 49% PBS, 50% Glycerol
Purity	Caprylic Acid - Ammonium Sulfate precipitation
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab193620** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/50 - 1/100. Predicted molecular weight: 30 kDa.

Target

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors