

## Product datasheet

# HRP Anti-Myc tag antibody ab1326

★★★★★ 3 Abreviews 19 References

### Overview

---

<b>Product name</b>	HRP Anti-Myc tag antibody
<b>Description</b>	HRP Rabbit polyclonal to Myc tag
<b>Host species</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Specificity</b>	In a direct ELISA against 1 ug of C-Myc peptide coupled to BSA, conjugate diluted to 1/20,000 produced a 1.0 OD in a 15 minute reaction with TMB as substrate.
<b>Tested applications</b>	<b>Suitable for:</b> ELISA, WB, ICC
<b>Immunogen</b>	Synthetic peptide corresponding to Myc tag conjugated to keyhole limpet haemocyanin. (Peptide available as <a href="#">ab13837</a> )

 [Run BLAST with](#)

 [Run BLAST with](#)

### General notes

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 7.00 Preservative: 0.002% Thimerosal (merthiolate) Constituents: 1.19% HEPES, 0.58% Sodium chloride, 0.2% BSA
<b>Purity</b>	Affinity purified
<b>Purification notes</b>	Antibodies were immunoaffinity purified using the peptide conjugated to a solid phase support, and conjugated to horseradish peroxidase.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

## Applications

---

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab1326 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/5000 - 1/20000.
WB	★★★★★ (3)	1/1000 - 1/5000. Can be blocked with <a href="#">Human c-Myc peptide (ab13837)</a> .
ICC		1/200 - 1/500.

## Target

---

**Relevance** Epitope tags are short peptide sequences that are easily recognized by tag-specific antibodies. Due to their small size, epitope tags do not affect the tagged protein's biochemical properties. Most often sequences encoding the epitope tag are included with target DNA at the time of cloning to produce fusion proteins containing the epitope tag sequence. This allows anti-epitope tag antibodies to serve as universal detection reagents for any tag containing protein produced by recombinant means. This means that anti-epitope tag antibodies are a useful alternative to generating specific antibodies to identify, immunoprecipitate or immunoaffinity purify a recombinant protein. The anti-epitope tag antibody is usually functional in a variety of antibody-dependent experimental procedures. Expression vectors producing epitope tag fusion proteins are available for a variety of host expression systems including bacteria, yeast, insect and mammalian cells.

**Cellular localization** Nuclear

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors