

Product datasheet

Human coronavirus SARS-CoV-2 spike glycoprotein peptide ab273063

Description

Product name	Human coronavirus SARS-CoV-2 spike glycoprotein peptide
Accession	QHD43416
Animal free	No
Nature	Synthetic
Species	Human coronavirus
Additional sequence information	(C terminal). 20 amino acid peptide within the last 50 amino acids of the SARS CoV-2 spike glycoprotein

Specifications

Our [Abpromise guarantee](#) covers the use of **ab273063** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Blocking - Blocking peptide for Anti-SARS-CoV-2 spike glycoprotein antibody - Coronavirus (ab272504)
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
------------------------------	--

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors