

Product datasheet

Human CPNE1 knockout HEK293T cell lysate ab258829

1 Image

Overview

<b>Product name</b>	Human CPNE1 knockout HEK293T cell lysate
<b>Product overview</b>	Knockout cell lysate achieved by CRISPR/Cas9.
<b>Parental Cell Line</b>	HEK293T
<b>Organism</b>	Human
<b>Mutation description</b>	Knockout achieved by using CRISPR/Cas9, Homozygous: 26 bp deletion in exon 2.
<b>Passage number</b>	<20
<b>Knockout validation</b>	Sanger Sequencing
<b>Reconstitution notes</b>	To use as WB control, resuspend the lyophilizate in 50 µL of LDS* Sample Buffer to have a final concentration of 2 mg/ml. For reducing conditions, we recommend a final concentration of 0.1 M DTT.

*\*Usage of SDS sample buffer is not recommended with these lyophilized lysates.*

**Notes**

**Lysate preparation:** Our lysates are made using RIPA buffer to which we add a protease inhibitor cocktail and phosphatase inhibitor cocktail (ratio: 300:100:10). *This means that the protein of interest is denatured.* If you require a native form of the protein please use the live cell version - found [here](#). Please refer to our lysis protocol for further details on how our lysates are prepared.

**User storage instructions:** After reconstitution, store the lysate at -80°C.

Access thousands of knockout cell lysates, generated from commonly used cancer cell lines. [See here for more information on knockout cell lysates.](#)

Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances. It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.

This product is subject to limited use licenses from The Broad Institute and ERS Genomics Limited, and is developed with patented technology. For full details of the limited use licenses and relevant patents please refer to our [limited use license](#) and [patent pages](#).

## Properties

**Storage instructions** Store at -80°C. Please refer to protocols.

Components	1 kit
ab261361 - Human CPNE1 knockout HEK293T cell lysate (Lyophilized)	1 x 100µg
ab255553 - Human wild-type HEK293T cell lysate (Lyophilized)	1 x 100µg

**Cell type** epithelial

**STR Analysis** Amelogenin X D5S818: 8, 9 D13S317: 12, 14 D7S820: 11 D16S539: 9, 13 vWA: 16, 19 TH01: 7, 9.3 TPOX: 11 CSF1PO: 11, 12

## Target

**Function** May function in membrane trafficking. Exhibits calcium-dependent phospholipid binding properties.

**Sequence similarities** Belongs to the copine family.  
Contains 2 C2 domains.  
Contains 1 VWFA domain.

## Images

Mut	CCCTCAGATGGCCCACTGGGTGACCTTG-----TTCAT
WT	CCCTCAGATGGCCCACTGGGTGACCTTGGTTCAGCTGTCCATTCCTGTGACCATCTCAT

Sanger Sequencing - Human CPNE1 knockout  
HEK293T cell lysate (ab258829)

Homozygous: 26 bp deletion in exon 2

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors