

Product datasheet

Human p53 DINP1/TP53INP1 peptide ab9776

1 References

Description

Product name	Human p53 DINP1/TP53INP1 peptide
Animal free	No
Nature	Synthetic
Species	Human
Sequence	MFQRLNKMFGVEVS amino acids 1 to 14 of human p53DINP1 (1).
Amino acids	1 to 14

Specifications

Our [Abpromise guarantee](#) covers the use of **ab9776** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Blocking - Blocking peptide for Anti-p53 DINP1/TP53INP1 antibody (ab9775)
Form	Liquid
Additional notes	This product was previously labelled as p53 DINP1

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle. pH: 7.20 Preservative: 0.02% Sodium azide Constituent: 0.1% BSA
------------------------------	---

General Info

Function	In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis.
Tissue specificity	Ubiquitously expressed.
Cellular localization	Nucleus > PML body.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors