

Product datasheet

Human Pallidin peptide ab23186

Overview

Product name Human Pallidin peptide

Description

Nature Synthetic

Amino Acid Sequence

Species Human

Sequence C-EREKQLTARPAKRM

Amino acids 159 to 172

Specifications

Our [Abpromise guarantee](#) covers the use of **ab23186** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Blocking - Blocking peptide for Anti-Pallidin antibody ([ab10145](#))

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

General Info

Function Involved in the development of lysosome-related organelles, such as melanosomes and platelet-dense granules. May play a role in intracellular vesicle trafficking, particularly in the vesicle-docking and fusion process.

Tissue specificity Widely expressed.

Post-translational modifications Phosphorylated.

Cellular localization Cytoplasm. Endomembrane system. It can exist as a soluble protein as well as a peripheral membrane protein.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors