

## Product datasheet

# Human PAMCI peptide ab45455

### Overview

---

**Product name** Human PAMCI peptide

### Description

---

**Nature** Synthetic

**Amino Acid Sequence**

**Species** Human

**Sequence** C-QDSETTVGDVLLST

**Amino acids** 421 to 435

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab45455** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Blocking

**Form** Liquid

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

### General Info

---

**Function** May play a role in regulating vesicular trafficking in cells.

**Sequence similarities** Contains 1 Ras-associating domain.

**Cellular localization** Endosome. Accumulates on perinuclear endosomes.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors