

Product datasheet

Human Resistin peptide ab42348

Overview

Product name Human Resistin peptide

Description

Nature Synthetic

Amino Acid Sequence

Species Human

Amino acids 33 to 90

Specifications

Our [Abpromise guarantee](#) covers the use of **ab42348** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Blocking

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles.

Constituent: 0.1% Acetic acid

General Info

Function Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.

Tissue specificity Expressed only in fatty tissues.

Sequence similarities Belongs to the resistin/FIZZ family.

Cellular localization Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors