

Product datasheet

Human TAB1 peptide ab39870

Overview

Product name Human TAB1 peptide

Description

Nature Synthetic

Amino Acid Sequence

Species Human

Specifications

Our [Abpromise guarantee](#) covers the use of **ab39870** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Biological activity The ab39870 peptide blocks antibodies that recognise 13 amino acids near the centre of human TAB1.

Applications Western blot
Blocking - Blocking peptide for Anti-TAB1 antibody ([ab25878](#))

Form Liquid

Additional notes The ab39870 peptide blocks antibodies that recognise 13 amino acids near the centre of human TAB1.

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Preservative: 0.02% Sodium Azide
Constituents: 0.1% BSA, PBS, pH 7.2

General Info

Function May be an important signaling intermediate between TGFB receptors and MAP3K7/TAK1. May play an important role in mammalian embryogenesis.

Tissue specificity Ubiquitous.

Sequence similarities Contains 1 PP2C-like domain.

Post-translational modifications

Monoubiquitinated. Deubiquitinated by Y.enterocolitica YopP.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors