

## Product datasheet

# Anti-IgG antibody [A set of murine monoclonal antibodies] ab210935

### 1 References

#### Overview

<b>Product name</b>	Anti-IgG antibody [A set of murine monoclonal antibodies]
<b>Description</b>	Mouse monoclonal [A set of murine monoclonal antibodies] to IgG
<b>Host species</b>	Mouse
<b>Tested applications</b>	<b>Suitable for:</b> ELISA, WB
<b>Species reactivity</b>	<b>Reacts with:</b> Rat <b>Does not react with:</b> Mouse, Goat, Human
<b>Immunogen</b>	Full length protein corresponding to Rat IgG.

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: 0.09% Sodium azide Constituent: 99% PBS
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	A set of murine monoclonal antibodies
<b>Isotype</b>	IgG

#### Applications

Our [Abpromise guarantee](#) covers the use of **ab210935** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/1000 - 1/5000.

Application	Abreviews	Notes
WB		1/500 - 1/10000. Predicted molecular weight: 36 kDa.

---

## Target

**Cellular localization**                      Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors