

Product datasheet

Anti-IL-3 antibody [EPR19870-40] ab222650

Recombinant RabMAb

Overview

| | |
|----------------------------|---|
| Product name | Anti-IL-3 antibody [EPR19870-40] |
| Description | Rabbit monoclonal [EPR19870-40] to IL-3 |
| Host species | Rabbit |
| Tested applications | Suitable for: Sandwich ELISA |
| Species reactivity | Reacts with: Mouse |
| Immunogen | The details of the immunogen for this antibody are not available. |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| Storage buffer | Preservative: 0.02% Sodium azide Constituent: PBS |
| Purity | Protein A purified |
| Clonality | Monoclonal |
| Clone number | EPR19870-40 |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab222650** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|----------------|-----------|--|
| Sandwich ELISA | | Use at an assay dependent concentration. |

Target

| | |
|------------------------------|--|
| Function | Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. This CSF induces granulocytes, macrophages, mast cells, stem cells, erythroid cells, eosinophils and megakaryocytes. |
| Tissue specificity | Activated T-cells, mast cells, natural killer cells. |
| Sequence similarities | Belongs to the IL-3 family. |
| Cellular localization | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors