



Product datasheet

Anti-IL-36R antibody [ABM47A2] - N-terminal ab210933

1 Image

Overview

Product name	Anti-IL-36R antibody [ABM47A2] - N-terminal
Description	Mouse monoclonal [ABM47A2] to IL-36R - N-terminal
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment corresponding to Human IL-36R aa 1-203 (N terminal). Sequence: MWSLLLCGLSIALPLSVTADGCKDIFMKNEILSASQPF AFNCTFPFITSG EVSVTWYKNSSKIPVSKIQRSRIHQDETWILFLPMEWGD SGVYQCVIKGR DSCHRIHVNLTVFEKHWCDTSIGGLPNLSDEYKQILHL GKDDSLTCHLHF PKSCVLGPIKWKDCNEIKGERFTVLETRLLVSNVSAE DRGNYACQAILT HSG Database link: Q9HB29
	 Run BLAST with  Run BLAST with
Positive control	Human small intestine lysate.
General notes	This product was previously labelled as IL-1RL2

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.05% Sodium azide Constituents: 99% PBS, 0.05% BSA
Purity	Protein G purified
Clonality	Monoclonal

Clone number	ABM47A2
Isotype	IgG2b
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab210933** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 2 µg/ml. Predicted molecular weight: 65 kDa.

Target

Function	Receptor for interleukin 1 family member 9 (IL1F9). Binding to the agonist leads to the activation of NF-kappa-B.
Sequence similarities	Belongs to the interleukin-1 receptor family. Contains 3 Ig-like C2-type (immunoglobulin-like) domains. Contains 1 TIR domain.
Cellular localization	Membrane.

Images



Anti-IL-36R antibody [ABM47A2] - N-terminal (ab210933) at 2 µg/ml + human small intestine lysate

Predicted band size: 65 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Apromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors