

Product datasheet

Anti-ILT-4 antibody [EPR22081] ab224701

Recombinant RabMAb

2 Images

Overview

Product name	Anti-ILT-4 antibody [EPR22081]
Description	Rabbit monoclonal [EPR22081] to ILT-4
Host species	Rabbit
Tested applications	Suitable for: ICC/IF, Flow Cyt Unsuitable for: IHC-P or WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment within Human ILT-4 aa 1-250. The exact sequence is proprietary. Database link: Q8N423
Positive control	ICC/IF: Human PBMC. Flow Cyt: Human PBMC.
General notes	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol, 0.05% BSA
Purity	Protein A purified
Clonality	Monoclonal
Clone number	EPR22081
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab224701** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
ICC/IF		1/50.
Flow Cyt		1/500.

Application notes Is unsuitable for IHC-P or WB.

Target

Function	Receptor for class I MHC antigens. Recognizes a broad spectrum of HLA-A, HLA-B, HLA-C and HLA-G alleles. Involved in the down-regulation of the immune response and the development of tolerance. Competes with CD8A for binding to class I MHC antigens. Inhibits FCGR1A-mediated phosphorylation of cellular proteins and mobilization of intracellular calcium ions.
Tissue specificity	Expressed on monocytes and B-cells, and at lower levels on dendritic cells. Detected at low levels in natural killer (NK) cells.
Sequence similarities	Contains 4 Ig-like C2-type (immunoglobulin-like) domains.
Domain	Contains 3 copies of a cytoplasmic motif that is referred to as the immunoreceptor tyrosine-based inhibitor motif (ITIM). This motif is involved in modulation of cellular responses. The phosphorylated ITIM motif can bind the SH2 domain of several SH2-containing phosphatases.
Post-translational modifications	Phosphorylated on tyrosine residues. Dephosphorylated by PTPN6.
Cellular localization	Membrane.

Images




Immunocytochemistry/ Immunofluorescence - Anti-ILT-4 antibody [EPR22081] (ab224701)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized Human PBMC (human primary peripheral blood mononuclear cell) cells labeling ILT-4 with ab224701 at 1/50 dilution, followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) (ab150077) secondary antibody at 1/1000 dilution (green). Confocal image showing co-staining with CD14 and cytoplasmic staining in PBMC cells. The nuclear counterstain is DAPI (blue). Counterstained with an Alexa Fluor® 647 anti-human CD14 antibody at a 1/100 dilution (red).

The **negative control** is the secondary antibody only.

ILT-4 is expressed in monocytes (PMID:10879687) and CD14 is an established marker for monocytes and macrophages (PMID: 23382732).



Flow Cytometry - Anti-ILT-4 antibody [EPR22081] (ab224701)

Flow cytometric analysis of human primary peripheral blood mononuclear cell (PBMC) labeling ILT-4 with ab224701 at 1/500 (right panel) compared with a Rabbit IgG, monoclonal [EPR25A] - Isotype Control (ab172730) (left panel). Goat anti rabbit IgG (Dylight® 488, ab98462) at 1/2000 dilution was used as the secondary antibody.

Cells were stained with Alexa Fluor® 647-conjugated CD14 and rabbit IgG (Left) or ILT-4 (Right).

Data shown were gated on viable cells.

The CD14 co-staining result observed is consistent with what has been described in the literature (PMID:10879687).

ILT-4 is expressed in monocytes and CD14 is an established marker for monocytes and macrophages (PMID: 23382732).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors