

## Product datasheet

# Immunohistochemistry - Methods Express Book Series ab94964

★★★★★ 1 Abreviews 1 Image

### Overview

---

**Product name** Immunohistochemistry - Methods Express Book Series

**General notes** **Edited and contributed to by four Abcam employees**

Immunohistochemistry has long been a cornerstone technique in biomedical research, but there have been major new developments in recent years, including a growing trend towards automation and high throughput technology. These new techniques, often coupled with some of the established methodology, have seen an expanding range of applications.

Immunohistochemistry: Methods Express is a comprehensive research guide that describes both the key new techniques and more established methods. Every chapter discusses the merits and limitations of the various approaches and then provides selected tried-and-tested protocols, as well as a plethora of good practical advice, for immediate use at the bench. This book is an essential laboratory manual for all life science researchers who use immunochemistry. It is a book which both experienced researchers and those new to the field will find invaluable.

**Contents:**

1. Introduction
2. Antibodies for immunochemistry
3. The selection of reporter labels
4. Immunochemical staining technique
5. Multiple immunochemical staining techniques
6. Confocal microscopy and immunohistochemistry
7. Ultrastructural immunochemistry
8. Image capture, analysis and quantification
9. Quality assurance in immunohistochemistry
10. Automated immunochemistry

Paperback: 300 pages

Publisher: Scion Publishing Ltd (15 Dec 2006)

Language English

Free shipping


### Properties

---

**Form** No information available

## Images

---



**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors