

Product datasheet

Anti-intestinal alkaline phosphatase antibody [AP1B9] ab18110

1 References

Overview

Product name	Anti-intestinal alkaline phosphatase antibody [AP1B9]
Description	Mouse monoclonal [AP1B9] to intestinal alkaline phosphatase
Host species	Mouse
Tested applications	Suitable for: IHC-P, IHC-Fr
Species reactivity	Reacts with: Cow
Immunogen	Full length native protein (purified) corresponding to Cow intestinal alkaline phosphatase.
General notes	Can be used as a second layer APAAP reagent for immunohistology.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Constituent: PBS
Purity	Protein A purified
Primary antibody notes	Can be used as a second layer APAAP reagent for immunohistology.
Clonality	Monoclonal
Clone number	AP1B9
Isotype	IgG1

Applications

Our [Abpromise guarantee](#) covers the use of **ab18110** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
-------------	-----------	-------

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

Target

Sequence similarities	Belongs to the alkaline phosphatase family.
Cellular localization	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors