

Product datasheet

Jurkat nuclear extract lysate (TPA + CI stimulated) ab14850

Overview

Product name	Jurkat nuclear extract lysate (TPA + CI stimulated)
General notes	<p>Cell line: Jurkat (Human T cell lymphoblast-like).</p> <p>Extracts have been quality control tested by Western blot and the Electrophoretic Mobility Shift Assay (EMSA).</p> <p>Cells were cultured in medium supplemented with 50 ng/mL TPA (phorbol, 12-myristate, 13 acetate) and 0.5 μM calcium ionophore A23187 (CI) for two hours at 37°C immediately prior to harvesting.</p>

Properties

Mycoplasma free	Yes
Form	Liquid
Storage instructions	Aliquot and store at -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	<p>pH: 7.90</p> <p>Constituents: 0.75% Potassium chloride, 0.009% PMSF, 0.01% Magnesium chloride, 0.008% DTT, 0.48% HEPES, 20% Glycerol</p>
Background	<p>Often this cell line is called 'JM' (Jurkat and JM are derived from the same patient and are sister clones), occasionally JM may be a subclone with somewhat divergent features. Jurkat cells can be transfected and are therefore useful for studies of blood proto-oncogenes expression, apoptosis and cell survival (e.g HIV). They produce IL-2 and can also be used for studies of differentiation.</p>

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors