

Product datasheet

Anti-Kallikrein 5 antibody ab7283

4 References

Overview

Product name	Anti-Kallikrein 5 antibody
Description	Rabbit polyclonal to Kallikrein 5
Host species	Rabbit
Immunogen	Synthetic peptide, corresponding to amino acids 280-293 of Human KLK5.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	<p>pH: 7.4</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p>
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Target

Function	May be involved in desquamation.
Tissue specificity	Expressed in skin, breast, brain and testis. Expressed at the stratum granulosum of palmar skin.
Sequence similarities	<p>Belongs to the peptidase S1 family. Kallikrein subfamily.</p> <p>Contains 1 peptidase S1 domain.</p>

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors