



Product datasheet

Anti-Keratinocyte differentiation associated protein antibody ab204583

[4 Images](#)

Overview

Product name	Anti-Keratinocyte differentiation associated protein antibody
Description	Rabbit polyclonal to Keratinocyte differentiation associated protein
Host species	Rabbit
Tested applications	Suitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	<p>Recombinant fragment corresponding to Human Keratinocyte differentiation associated protein aa 24-81.</p> <p>Sequence:</p> <p>TLGGPEEEESTIENYASRPEAFNTPFLNIDKLRSFAKAEFL NWHALFESI KRKLPFLN</p> <p>Database link: P60985</p> <p style="text-align: right;">  Run BLAST with  Run BLAST with </p>
Positive control	IHC-P: Human cervix, skin tissue.
General notes	<p>Reproducibility is key to advancing scientific discovery and accelerating scientists' next breakthrough.</p> <p>Abcam is leading the way with our range of recombinant antibodies, knockout-validated antibodies and knockout cell lines, all of which support improved reproducibility.</p> <p>We are also planning to innovate the way in which we present recommended applications and species on our product datasheets, so that only applications & species that have been tested in our own labs, our suppliers or by selected trusted collaborators are covered by our Abpromise™ guarantee.</p> <p>In preparation for this, we have started to update the applications & species that this product is Abpromise guaranteed for.</p> <p>We are also updating the applications & species that this product has been “predicted to work with,” however this information is not covered by our Abpromise guarantee.</p> <p>Applications & species from publications and Abreviews that have not been tested in our own labs or in those of our suppliers are not covered by the Abpromise guarantee.</p> <p>Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team</p>

ahead of purchase. Recommended alternatives for this product can be found below, as well as customer reviews and Q&As.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol, 59% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab204583** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/1000 - 1/2500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Target

Function	May act as a soluble regulator of keratinocyte differentiation. May play an important role in embryonic skin morphogenesis.
Tissue specificity	Highly expressed in skin and detected at lower levels in thymus. In skin, found exclusively in lamellar granules of granular keratinocytes and in the intracellular space of the stratum corneum. Also highly expressed in oral mucosa, tongue, esophagus, and stomach, and at much lower levels in bladder and uterus. Not detected in gastrointestinal mucosa.
Cellular localization	Secreted.


Images



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Keratinocyte differentiation associated protein antibody (ab204583)

Paraffin embedded human cervix tissue stained for Keratinocyte differentiation associated protein using ab204583 at 1/2500 dilution in immunohistochemical analysis.


Heat mediated antigen retrieval was performed with citrate buffer pH 6 before the IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Keratinocyte differentiation associated protein antibody (ab204583)

Paraffin embedded human skin tissue stained for Keratinocyte differentiation associated protein using ab204583 at 1/2500 dilution in immunohistochemical analysis.


Heat mediated antigen retrieval was performed with citrate buffer pH 6 before the IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Keratinocyte differentiation associated protein antibody (ab204583)

Paraffin embedded human pancreas tissue stained for Keratinocyte differentiation associated protein using ab204583 at 1/2500 dilution in immunohistochemical analysis. Shows no positivity in exocrine glandular cells as expected.

Heat mediated antigen retrieval was performed with citrate buffer pH 6 before the IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Keratinocyte differentiation associated protein antibody (ab204583)

Paraffin embedded human liver tissue stained for Keratinocyte differentiation associated protein using ab204583 at 1/2500 dilution in immunohistochemical analysis. Shows no positivity in hepatocytes cells as expected.

Heat mediated antigen retrieval was performed with citrate buffer pH 6 before the IHC staining protocol.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors