

## Product datasheet

# Anti-KRAP antibody ab195334

1 Image

### Overview

<b>Product name</b>	Anti-KRAP antibody
<b>Description</b>	Rabbit polyclonal to KRAP
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> IP
<b>Species reactivity</b>	<p><b>Reacts with:</b> Human</p> <p><b>Predicted to work with:</b> Sheep, Rabbit, Cat, Dog, Chimpanzee, Macaque monkey, Rhesus monkey, Gorilla, Common marmoset, Orangutan </p>
<b>Immunogen</b>	<p>Synthetic peptide within Human KRAP aa 1209-1259. The exact sequence is proprietary. NP_001123917.1.</p> <p>Sequence:</p> <p>MDRPLSSSAEAELEELEWQVASRRRKAWAKCRSSWQAS ETEDLSTEATTQD</p> <p>Database link: <a href="#">P28290</a></p> <p style="text-align: right;"> <a href="#">Run BLAST with</a>       <a href="#">Run BLAST with</a> </p>
<b>Positive control</b>	HeLa and 293T whole cell lysates.
<b>General notes</b>	<p>Previously labelled as SSFA2.</p> <p>Reproducibility is key to advancing scientific discovery and accelerating scientists' next breakthrough.</p> <p>Abcam is leading the way with our range of recombinant antibodies, knockout-validated antibodies and knockout cell lines, all of which support improved reproducibility.</p> <p>We are also planning to innovate the way in which we present recommended applications and species on our product datasheets, so that only applications &amp; species that have been tested in our own labs, our suppliers or by selected trusted collaborators are covered by our Abpromise™ guarantee.</p> <p>In preparation for this, we have started to update the applications &amp; species that this product is Abpromise guaranteed for.</p> <p>We are also updating the applications &amp; species that this product has been “predicted to work with,” however this information is not covered by our Abpromise guarantee.</p> <p>Applications &amp; species from publications and Abreviews that have not been tested in our own labs or in those of our suppliers are not covered by the Abpromise guarantee.</p> <p>Please check that this product meets your needs before purchasing. If you have any questions,</p>

special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, as well as customer reviews and Q&As.

## Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7 Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate
<b>Purity</b>	pH 7-8 Immunogen affinity purified
<b>Purification notes</b>	ab195334 was affinity purified using an epitope specific to KRAP immobilized on a solid support.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

## Applications

---

Our [Abpromise guarantee](#) covers the use of **ab195334** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IP		Use at 2-10 µg/mg of lysate.


---

## Target

---

## Images

---



Immunoprecipitation of 293T whole cell lysate in NETN lysis buffer using ab195334 (lane 1) at 6  $\mu\text{g}/\text{mg}$  lysate (0.5 or 1 mg/IP; 20% of IP loaded/lane). Immunoblot analysis of KRAP utilized ab195334 at 1  $\mu\text{g}/\text{mL}$ . The lysate in lane 2 is precipitated with an IgG control antibody. Detection utilized chemiluminescence with an exposure time of 3 minutes.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors