

Product datasheet

Anti-Malondialdehyde antibody ab27644

2 References

Overview

Product name	Anti-Malondialdehyde antibody
Description	Goat polyclonal to Malondialdehyde
Host species	Goat
Specificity	This antibody specifically binds to Malondialdehyde LDL and other MDA modified proteins.
Immunogen	Chemical/ Small Molecule MDA modified protein
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
Storage buffer	<p>pH: 7.20</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 0.2% PBS, 0.0292% EDTA, 0.435% Sodium chloride, 30% Glycerol (glycerin, glycerine)</p>
Purification notes	This antibody was purified by Malondialdehyde modified protein Sepharose affinity column.
Clonality	Polyclonal
Isotype	IgG

Target

Relevance	Malondialdehyde (MDA) is a natural product formed in all mammalian cells as a product of lipid peroxidation. MDA is a highly reactive three carbon dialdehyde produced as a byproduct of polyunsaturated fatty acid peroxidation and arachidonic acid metabolism. MDA readily combines
------------------	--

with several functional groups on molecules including proteins, lipoproteins, and DNA. It reacts with DNA to form adducts to deoxyguanosine and deoxyadenosine. The major adduct to DNA is a pyrimidopurinone called M1G which appears to be a major endogenous DNA adduct in human beings that may contribute significantly to cancer linked to lifestyle and dietary factors. MDA modified proteins may show altered physico chemical behavior and antigenicity. MDA is toxic and has been implicated in aging mutagenesis, carcinogenesis, diabetic nephropathy and radiation damage. Increased expression of MDA has been reported in the brains of Alzheimer's patients. Antibodies to MDA will help to visualize the MDA adducts.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors