

Product datasheet

Anti-MAML3 antibody [797C2a] ab53221

1 Image

Overview

Product name	Anti-MAML3 antibody [797C2a]
Description	Mouse monoclonal [797C2a] to MAML3
Host species	Mouse
Tested applications	Suitable for: WB, Dot blot
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment of human MAML3.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 1% BSA, PBS
Purity	Protein G purified
Purification notes	ab53221 was purified from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum and filtered through a 0.22 µm membrane.
Clonality	Monoclonal
Clone number	797C2a
Isotype	IgG1

Applications

Our [Abpromise guarantee](#) covers the use of **ab53221** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes

Application	Abreviews	Notes
WB		Use at an assay dependent dilution. Detects a band of approximately 34 kDa (predicted molecular weight: 122 kDa).
Dot blot		Use at an assay dependent dilution.

Target

Images



Anti-MAML3 antibody [797C2a] (ab53221) + immunised recombinant protein fragment

Predicted band size: 122 kDa

Observed band size: 34 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors