

Product datasheet

MES Run Buffer (20X) - RunBlue™ (Bis-Tris) ab270224

1 Image

Overview

Product name

MES Run Buffer (20X) - RunBlue™ (Bis-Tris)

General notes

MES Run Buffer (20X) - RunBlue™ (Bis-Tris) (ab270224) is designed for use with Bis-Tris Protein Gels. This buffer is made using high purity ingredients. Buffers are provided at 20X concentration, and it is recommended to dilute with deionized water for use.

This product is manufactured by Expedeon, an Abcam company. Expedeon product code NXB70425 was previously called RunBlue MES Run Buffer - 4L (20x Concentrate) and is the same as the 4000 mL size of this product. Expedeon product code NXB70500 was previously called RunBlue MES Run Buffer - 500 ml (20x Concentrate) and is the same as the 500 mL size.

ab270224 is recommended for small to medium sized proteins run on any RunBlue™ Bis-Tris gel. Using pre-mixed buffers is a convenient way to decrease variability between results and save time in the lab.

Please note that for TEO-Tricine Precast Gels – RunBlue™ we recommend using Run Buffer (20X) -RunBlue™ (TEO-Tricine) ([ab119197](#)).

Properties

Form

Liquid

Storage instructions


Shipped at Room Temperature. Store at 4°C or -20°C.

Storage buffer

pH: 7.3

Constituent: 5% Sodium lauryl sulfate

Images



Protein migration chart.

SDS-PAGE - MES Run Buffer (20X) - RunBlue™
(Bis-Tris) (ab270224)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors