

Product datasheet

Anti-Mipu1 antibody ab106432

3 Images

Overview

Product name	Anti-Mipu1 antibody
Description	Rabbit polyclonal to Mipu1
Host species	Rabbit
Tested applications	Suitable for: ICC/IF, IHC-P, WB
Species reactivity	Reacts with: Mouse, Rat, Human
Immunogen	Synthetic peptide corresponding to Human Mipu1 (N terminal). Database link: Q5HYK9
Positive control	WB: Human heart tissue lysate IHC-P: Human brain tissue
General notes	This product was previously labelled as ZNF667

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.02% Sodium azide Constituent: PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab106432** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
ICC/IF		Use a concentration of 20 µg/ml.

Application	Abreviews	Notes
IHC-P		Use a concentration of 2.5 - 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 70 kDa.

Target


Function	May be involved in transcriptional regulation.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 15 C2H2-type zinc fingers. Contains 1 KRAB domain.
Cellular localization	Nucleus.

Images




Anti-Mipu1 antibody (ab106432) at 1 µg/ml + Human heart tissue lysate

Predicted band size: 70 kDa



ab106432, at 2.5 µg/ml, staining Mipu1 in paraffin embedded Human brain tissue by Immunohistochemistry.



Immunohistochemical analysis of human brain tissue labelling
Mipu1 with ab106432 at 20 ug/mL.

Immunocytochemistry/ Immunofluorescence - Anti-
Mipu1 antibody (ab106432)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors