

Product datasheet

MOPS Run Buffer (20X) - RunBlue™ (Bis-Tris) ab270225

1 Image


Overview

Product name	MOPS Run Buffer (20X) - RunBlue™ (Bis-Tris)
General notes	<p>MOPS Run Buffer (20X) - RunBlue™ (Bis-Tris) (ab270225) is designed for use with Bis-Tris Protein Gels. This buffer is made using high purity ingredients. Buffers are provided at 20X concentration, and it is recommended to dilute with deionized water for use.</p> <p>This product is manufactured by Expedeon, an Abcam company. Expedeon product code NXB75425 was previously called RunBlue MOPS Run Buffer - 4L (20x Concentrate) and is the same as the 4000 mL size of this product. Expedeon product code NXB75500 was previously called RunBlue MOPS Run Buffer - 500 ml (20x Concentrate) and is the same as the 500 mL size.</p> <p>ab270225 is recommended for medium to large sized proteins run on any RunBlue™ Bis-Tris gel. Using pre-mixed buffers is a convenient way to decrease variability between results and save time in the lab.</p> <p>Please note that for TEO-Tricine Precast Gels – RunBlue™ we recommend using Run Buffer (20X) -RunBlue™ (TEO-Tricine) (ab119197).</p>

Properties

Form	Liquid
Storage instructions	Shipped at Room Temperature. Store at Room Temperature.
Storage buffer	pH: 7.8 Constituent: 5% Sodium lauryl sulfate

Images



Protein migration chart.

SDS-PAGE - MOPS Run Buffer (20X) - RunBlue™
(Bis-Tris) (ab270225)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors