

Product datasheet

Mouse IgM, kappa monoclonal [MM-30] - isotype control ab18401

6 References

Overview

Product name Mouse IgM, kappa monoclonal [MM-30] - isotype control

Tested applications **Suitable for:** IHC-P, IHC-Fr, WB, IP, Flow Cyt

Immunogen KLH

General notes This immunoglobulin reacts with KLH.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

Properties

Form Liquid

Storage instructions Shipped at 4°C. Store at +4°C.

Storage buffer pH: 7.20
Preservative: 0.09% Sodium azide
Constituent: PBS

Purification notes The immunoglobulin was purified by affinity chromatography.

Isotype control notes This immunoglobulin reacts with KLH.

Clonality Monoclonal

Clone number MM-30

Isotype IgM

Light chain type kappa

Applications

The Abpromise guarantee

Our [Abpromise guarantee](#) covers the use of ab18401 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. Suitable for surface or intracellular staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors