

Product datasheet

Mouse IL-18 Antibody Pair - BSA and Azide free
ab241886

Recombinant RabMAb

5 Images

Overview

Product name	Mouse IL-18 Antibody Pair - BSA and Azide free
Assay type	ELISA set
Range	31.25 pg/ml - 2000 pg/ml
Species reactivity	Reacts with: Mouse
Product overview	<p>The Antibody Pair can be used to quantify Mouse IL-18. BSA and Azide free antibody pairs include unconjugated capture and detector antibodies suitable for sandwich ELISAs. The antibodies are provided at an approximate concentration of 1 mg/ml as measured by the protein A280 method. The recommended antibody orientation is based on internal optimization for ELISA-based assays. Antibody orientation is assay dependent and needs to be optimized for each assay type. Both capture and detector antibodies are rabbit monoclonal antibodies delivering consistent, specific, and sensitive results.</p> <p>For additional information on the performance of the antibody pair, see the equivalent SimpleStep ELISA® Kit (ab216165), which uses the same antibodies. However, due to differences in their formulation, this antibody pair cannot be used with the consumables provided with our SimpleStep ELISA Kits. Please note that the range provided for the pairs is only an estimation based on the performance of the related product using the same antibody pair. Performance of the antibody pair will depend on the specific characteristics of your assay. We guarantee the product works in sandwich ELISA, but we do not guarantee the sensitivity or dynamic range of the antibody pair in your assay.</p> <p>Download SDS here.</p>
Tested applications	Suitable for: Sandwich ELISA
Platform	Reagents

Properties

Storage instructions	Store at +4°C. Please refer to protocols.
Carrier free	Yes

Components	10 x 96 tests
Mouse IL-18 Capture Antibody (unconjugated)	1 x 100µg
Mouse IL-18 Detector Antibody (unconjugated)	1 x 100µg

Function Augments natural killer cell activity in spleen cells and stimulates interferon gamma production in T-helper type I cells.

Sequence similarities Belongs to the IL-1 family.

Cellular localization Secreted.


Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab241886 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Sandwich ELISA		Use at an assay dependent concentration.

Images




Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab216165](#)), which uses the same antibody pair. For additional information on the performance of pair and kit, refer to the corresponding kit datasheet. Due to differences in the formulation and format of the antibodies in this pair, they cannot be used as substitutes for the antibody components in our SimpleStep ELISA® Kits.

Sandwich ELISA - Mouse IL-18 Antibody Pair - BSA and Azide free (ab241886)

Standard Curve Measurements			
Concentration (pg/mL)	O.D 450 nm		Mean O.D
	1	2	
0	0.055	0.055	0.055
31.3	0.103	0.101	0.102
62.5	0.150	0.150	0.150
125	0.246	0.250	0.248
250	0.422	0.432	0.427
500	0.779	0.820	0.799
1,000	1.469	1.534	1.501
2,000	2.651	2.755	2.703

Sandwich ELISA - Mouse IL-18 Antibody Pair - BSA and Azide free (ab241886)

Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab216165](#)), which uses the same antibody pair. For additional information on the performance of pair and kit, refer to the corresponding kit datasheet. Due to differences in the formulation and format of the antibodies in this pair, they cannot be used as substitutes for the antibody components in our SimpleStep ELISA® Kits.



Sandwich ELISA - Mouse IL-18 Antibody Pair - BSA and Azide free (ab241886)

Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab216165](#)), which uses the same antibody pair. For additional information on the performance of pair and kit, refer to the corresponding kit datasheet. Due to differences in the formulation and format of the antibodies in this pair, they cannot be used as substitutes for the antibody components in our SimpleStep ELISA® Kits.

Standard Curve Measurements			
Concentration (pg/mL)	O.D 450 nm		Mean O.D
	1	2	
0	0.063	0.065	0.064
31.3	0.113	0.114	0.113
62.5	0.162	0.193	0.177
125	0.261	0.261	0.261
250	0.455	0.450	0.452
500	0.831	0.862	0.847
1,000	1.515	1.568	1.542
2000	2.743	2.703	2.723

Sandwich ELISA - Mouse IL-18 Antibody Pair - BSA and Azide free (ab241886)

Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab216165](#)), which uses the same antibody pair. For additional information on the performance of pair and kit, refer to the corresponding kit datasheet. Due to differences in the formulation and format of the antibodies in this pair, they cannot be used as substitutes for the antibody components in our SimpleStep ELISA® Kits.

Powered by recombinant antibodies



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Sandwich ELISA - Mouse IL-18 Antibody Pair - BSA and Azide free (ab241886)

To learn more about the advantages of recombinant antibodies see [here](#).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors