

Product datasheet

Mouse monoclonal [5A9] Anti-Human IgG Fc (HRP) ab7499

1 References

Overview

Product name	Mouse monoclonal [5A9] Anti-Human IgG Fc (HRP)
Host species	Mouse
Target species	Human
Tested applications	Suitable for: Sandwich ELISA, ELISA
Immunogen	These clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunised with IgG from human serum.
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.40 Preservative: 0.05% Proclin 300 Constituents: 0.8% Sodium chloride, 0.02% Potassium chloride, 0.024% Monobasic dihydrogen potassium phosphate, 0.18% Dibasic monohydrogen sodium phosphate
Purity	Protein A purified
Clonality	Monoclonal
Clone number	5A9
Isotype	IgG2a
General notes	Concentration varies from lot to lot and can be provided on request.

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab7499 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Sandwich ELISA		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors