

Product datasheet

Anti-Mycoplasma synoviae antibody ab34145

Overview

Product name	Anti-Mycoplasma synoviae antibody
Description	Chicken polyclonal to Mycoplasma synoviae
Host species	Chicken
Tested applications	Suitable for: ELISA, Neutralising
Species reactivity	Reacts with: Other species
Immunogen	Mycoplasma synoviae

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Constituent: Whole serum
Purity	Whole antiserum
Purification notes	Sterile filtered serum produced in SPF (Specific Pathogen Free) chicken.
Clonality	Polyclonal
Isotype	IgY

Applications

Our [Abpromise guarantee](#) covers the use of **ab34145** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.
Neutralising		Use at an assay dependent dilution.

Target

Relevance

Mycoplasma synoviae is a common pathogen of poultry, with worldwide economical importance, that causes respiratory infections and synovitis. Infection occurs mainly in chickens and turkeys but many other avian species can be affected by this pathogen (ducks, geese, etc).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors