

Product datasheet

Anti-NAGLU/NAG antibody [EPR20708] ab214671

Recombinant RabMAb

1 References 11 Images

Overview

<b>Product name</b>	Anti-NAGLU/NAG antibody [EPR20708]
<b>Description</b>	Rabbit monoclonal [EPR20708] to NAGLU/NAG
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB, ICC/IF, Flow Cyt, IHC-P
<b>Species reactivity</b>	<b>Reacts with:</b> Mouse, Rat, Human
<b>Immunogen</b>	<b>This product was produced with the following immunogens:</b> Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.  Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
<b>Positive control</b>	WB: Human placenta, liver and fetal kidney lysates; HepG2 whole cell lysate; Mouse colon and testis lysates; Rat colon, testis, liver and spleen lysates. IHC-P: Human kidney and liver tissues; Mouse and rat kidney tissues. ICC/IF: HepG2 and Hepa1-6 cells. Flow Cyt: Hepa1-6 and HepG2 cells.
<b>General notes</b>	This product was previously labelled as NAGLU

This product is a recombinant monoclonal antibody, which offers several advantages including:

- High batch-to-batch consistency and reproducibility
- Improved sensitivity and specificity
- Long-term security of supply
- Animal-free production

For more information [see here](#).

Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb<sup>®</sup> patents](#).

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.2

	Preservative: 0.01% Sodium azide Constituents: 0.05% BSA, 40% Glycerol, PBS
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	EPR20708
<b>Isotype</b>	IgG

## Applications


**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab214671 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
<b>WB</b>		1/1000. Detects a band of approximately 82 kDa (predicted molecular weight: 82 kDa). Additional bands observed in some rodent tissue lysates where NAGLU is not highly expressed.
<b>ICC/IF</b>		1/100.
<b>Flow Cyt</b>		1/50.
<b>IHC-P</b>		1/4000. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.

## Target

<b>Function</b>	Involved in the degradation of heparan sulfate.
<b>Tissue specificity</b>	Liver, ovary, peripheral blood leukocytes, testis, prostate, spleen, colon, lung, placenta and kidney.
<b>Involvement in disease</b>	Defects in NAGLU are the cause of mucopolysaccharidosis type 3B (MPS3B) [MIM:252920]; also known as Sanfilippo syndrome B. MPS3B is a form of mucopolysaccharidosis type 3, an autosomal recessive lysosomal storage disease due to impaired degradation of heparan sulfate. MPS3 is characterized by severe central nervous system degeneration, but only mild somatic disease. Onset of clinical features usually occurs between 2 and 6 years; severe neurologic degeneration occurs in most patients between 6 and 10 years of age, and death occurs typically during the second or third decade of life.
<b>Cellular localization</b>	Lysosome.

## Images



**Lanes 1-3 :** Anti-NAGLU/NAG antibody [EPR20708] (ab214671) at 1/1000 dilution

**Lane 4 :** Anti-NAGLU/NAG antibody [EPR20708] (ab214671) at 1/5000 dilution

**Lane 1 :** Human placenta lysate at 20 µg

**Lane 2 :** Human liver lysate at 20 µg

**Lane 3 :** Human fetal kidney lysate at 10 µg

**Lane 4 :** HepG2 (human liver hepatocellular carcinoma cell line) whole cell lysate at 10 µg

### Secondary

**Lanes 1-2 :** VeriBlot for IP Detection Reagent (HRP) (ab131366) at 1/4000 dilution

**Lanes 3-4 :** Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution


Developed using the ECL technique.

**Predicted band size:** 82 kDa

**Observed band size:** 82 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.

Exposure time: Lanes 1-2: 1 minute; Lanes 3-4: 30 seconds.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

Immunohistochemical analysis of paraffin-embedded human kidney tissue labeling NAGLU/NAG with ab214671 at 1/4000 dilution, followed by Goat Anti-Rabbit IgG H&L (HRP) Ready to use.

Granular cytoplasmic staining on human kidney tubules (PMID: 4291567, PMID: 24244710). Counter stained with hematoxylin.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (HRP) Ready to use.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.




Immunocytochemistry/ Immunofluorescence - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

Immunofluorescent analysis of 100% methanol-fixed HepG2 (human liver hepatocellular carcinoma cell line) cells labeling NAGLU/NAG with ab214671 at 1/100 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) (ab150077) secondary antibody at 1/1000 dilution (green). Confocal image showing cytoplasmic staining on HepG2 cell line.

The nuclear counter stain is DAPI (blue). Tubulin is detected with Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) (ab195889) (red) at 1/200 dilution.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) (ab150077) secondary antibody at 1/1000 dilution.



Western blot - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

**Lanes 1-4** : Anti-NAGLU/NAG antibody [EPR20708] (ab214671) at 1/2000 dilution

**Lanes 5-6** : Anti-NAGLU/NAG antibody [EPR20708] (ab214671) at 1/1000 dilution

**Lane 1** : Mouse colon lysate at 20 µg

**Lane 2** : Rat testis lysate at 20 µg

**Lane 3** : Mouse testis lysate at 20 µg

**Lane 4** : Rat colon lysate at 20 µg

**Lane 5** : Rat liver lysate at 10 µg

**Lane 6** : Rat spleen lysate at 10 µg

### Secondary

**All lanes** : Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at

1/100000 dilution


Developed using the ECL technique.

**Predicted band size:** 82 kDa

**Observed band size:** 82 kDa

Blocking/Dilution buffer: 5% NFDN/TBST.

Exposure time: Lanes 1-2,5-6: 3 minutes; Lanes 3-4: 30 seconds.




Immunohistochemical analysis of paraffin-embedded human liver tissue labeling NAGLU/NAG with ab214671 at 1/4000 dilution, followed by Goat Anti-Rabbit IgG H&L (HRP) Ready to use.

Granular cytoplasmic staining on human liver (PMID: 4291567, PMID: 8776591). Counter stained with hematoxylin.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (HRP) Ready to use.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)




Immunohistochemical analysis of paraffin-embedded mouse kidney tissue labeling NAGLU/NAG with ab214671 at 1/4000 dilution, followed by Goat Anti-Rabbit IgG H&L (HRP) Ready to use.

Granular cytoplasmic staining on mouse kidney tubules (PMID: 4291567, PMID: 24244710). Counter stained with hematoxylin.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (HRP) Ready to use.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

Immunohistochemical analysis of paraffin-embedded rat kidney tissue labeling NAGLU/NAG with ab214671 at 1/4000 dilution, followed by Goat Anti-Rabbit IgG H&L (HRP) Ready to use. Granular cytoplasmic staining on rat kidney tubules (PMID: 4291567, PMID: 24244710). Counter stained with hematoxylin.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (HRP) Ready to use.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.




Immunocytochemistry/ Immunofluorescence - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

Immunofluorescent analysis of 100% methanol-fixed Hepa1-6 (mouse hepatoma epithelial cell line) cells labeling NAGLU/NAG with ab214671 at 1/100 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) (ab150077) secondary antibody at 1/1000 dilution (green). Confocal image showing cytoplasmic staining on Hepa1-6 cell line.


The nuclear counter stain is DAPI (blue). Tubulin is detected with Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) (ab195889) (red) at 1/200 dilution.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) (ab150077) secondary antibody at 1/1000 dilution.



Flow Cytometry (Intracellular) - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

Flow cytometric analysis of 4% paraformaldehyde-fixed Hepa1-6 (mouse hepatoma epithelial cell line) cell line labeling NAGLU/NAG with ab214671 at 1/50 dilution (red) compared with a Rabbit IgG, monoclonal [EPR25A] - Isotype Control (ab172730) (black) and an unlabeled control (cells without incubation with primary antibody and secondary antibody) (blue). Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) at 1/2000 dilution was used as the secondary antibody.



Flow cytometric analysis of 4% paraformaldehyde-fixed HepG2 (human liver hepatocellular carcinoma cell line) cell line labeling NAGLU/NAG with ab214671 at 1/50 dilution (red) compared with a Rabbit IgG, monoclonal [EPR25A] - Isotype Control (ab172730) (black) and an unlabeled control (cells without incubation with primary antibody and secondary antibody) (blue). Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) at 1/2000 dilution was used as the secondary antibody.

Flow Cytometry (Intracellular) - Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

**Why choose a recombinant antibody?**

 <b>Research with confidence</b> Consistent and reproducible results	 <b>Long-term and scalable supply</b> Recombinant technology
 <b>Success from the first experiment</b> Confirmed specificity	 <b>Ethical standards compliant</b> Animal-free production

Anti-NAGLU/NAG antibody [EPR20708] (ab214671)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors