

Product datasheet

Native Cow AGE-BSA protein ab129535

Description

Product name	Native Cow AGE-BSA protein
Expression system	Native
Protein length	Full length protein
Animal free	No
Nature	Native
Species	Cow

Specifications

Our [Abpromise guarantee](#) covers the use of **ab129535** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Blocking - Blocking peptide for Anti-AGE antibody (ab23722)
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Constituent: PBS
------------------------------	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors