

Product datasheet

Native Dog IgG protein ab198651

1 Image

Description

Product name	Native Dog IgG protein
Purity	> 95 % SDS-PAGE.
Expression system	Native
Protein length	Full length protein
Animal free	No
Nature	Native
Species	Dog
Predicted molecular weight	150 kDa
Additional sequence information	Prepared from dog plasma.

Specifications

Our [Abpromise guarantee](#) covers the use of **ab198651** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized


Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -80°C. Avoid freeze / thaw cycle. pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.33% Sodium phosphate, 0.87% Sodium chloride
------------------------------	---

Reconstitution	We suggest the addition of de-ionized H2O to the original volume, followed by gentle swirling and/or vortexing to ensure adequate homogenization. If further dilution is required, we suggest using the lyophilization buffer.
-----------------------	--

General Info

Images



SDS PAGE analysis of ab198651:

Lane 1: 5 µg, under non reduced and not heated conditions;

Lane 2: 10 µg, under non reduced and not heated conditions;

Lane 3: 20 µg, under non reduced and not heated conditions.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors