

Product datasheet

Native human Thrombin protein (Active) ab62452

1 Image

Description

Product name	Native human Thrombin protein (Active)
Biological activity	This product can be used for cleavage of fusion proteins. The activity varies form lot to lot. It is typically 3800 U/mg.
Purity	> 95 % SDS-PAGE.
Expression system	Native
Accession	P00734
Protein length	Full length protein
Animal free	No
Nature	Native
Species	Human

Specifications

Our [Abpromise guarantee](#) covers the use of **ab62452** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Functional Studies
Form	Liquid
Additional notes	Concentration varies from lot to lot but is between 5mg/ml and 10mg/ml.

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. Constituent: 50% Glycerol This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------------------------	---

General Info

Function	Thrombin, which cleaves bonds after Arg and Lys, converts fibrinogen to fibrin and activates
-----------------	--

factors V, VII, VIII, XIII, and, in complex with thrombomodulin, protein C. Functions in blood homeostasis, inflammation and wound healing.

Tissue specificity

Expressed by the liver and secreted in plasma.

Involvement in disease

Factor II deficiency
Ischemic stroke
Thrombophilia due to thrombin defect
Pregnancy loss, recurrent, 2

Sequence similarities

Belongs to the peptidase S1 family.
Contains 1 Gla (gamma-carboxy-glutamate) domain.
Contains 2 kringle domains.
Contains 1 peptidase S1 domain.


Post-translational modifications

The gamma-carboxyglutamyl residues, which bind calcium ions, result from the carboxylation of glutamyl residues by a microsomal enzyme, the vitamin K-dependent carboxylase. The modified residues are necessary for the calcium-dependent interaction with a negatively charged phospholipid surface, which is essential for the conversion of prothrombin to thrombin.
N-glycosylated. N-glycan heterogeneity at Asn-121: Hex3HexNAc3 (minor), Hex4HexNAc3 (minor) and Hex5HexNAc4 (major). At Asn-143: Hex4HexNAc3 (minor) and Hex5HexNAc4 (major).

Cellular localization

Secreted, extracellular space.

Images



Novex 4-12% Bis-Tris
1 µg per lane
MOPS

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors