

## Product datasheet

# Nitrocellulose Transfer Membrane- 0.22 µm ab133413

## 1 References

### Overview

---

<b>Product name</b>	Nitrocellulose Transfer Membrane- 0.22 µm
<b>General notes</b>	High quality 0.22 um nitrocellulose membranes that are more durable than standard nitrocellulose membranes. ab133413 yields clearer blots as the membrane is less prone to distortion. The membranes have been pre-cut for ease of use. Compatible with gels 10 x 10 cm and smaller.

### Properties

---

<b>Form</b>	Solid
<b>Storage instructions</b>	Store at room temperature.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors