



Product datasheet

Anti-NLRC3 antibody ab222200

1 Image

Overview

Product name	Anti-NLRC3 antibody
Description	Rabbit polyclonal to NLRC3
Host species	Rabbit
Tested applications	Suitable for: ICC/IF
Species reactivity	Reacts with: Human
Immunogen	<p>Recombinant fragment corresponding to Human NLRC3 aa 58-146. Sequence:</p> <p>SNDSRIQRHRKALLSKVGGGPELGGPWHRLASLLLVE GLTDLQLREHDFI QVEATRGGGHPARTVALDRLFLPLSRVSVPPRVSTITG V</p> <p>Database link: Q7RTR2</p> <p style="text-align: right;">  Run BLAST with  Run BLAST with </p>
Positive control	ICC/IF: NB-4 cells.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	<p>pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol, PBS</p>
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab222200** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use a concentration of 1 - 4 µg/ml. Fixation/Permeabilization: PFA/Triton X-100.

Target

Function	May modulate T-cell activation. Decreases the transcription of genes that are normally up-regulated after T-cell stimulation. Delays degradation of NFKBIA/IKBA.
Tissue specificity	Detected in peripheral blood mononuclear cells.
Sequence similarities	Belongs to the NLRP family. Contains 16 LRR (leucine-rich) repeats. Contains 1 NACHT domain.
Cellular localization	Cytoplasm.

Images



PFA-fixed, Triton X-100 permeabilized NB-4 cells stained for NLRC3 (green) using ab222200 at 4 µg/ml in ICC/IF.

Immunocytochemistry/ Immunofluorescence - Anti-NLRC3 antibody (ab222200)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors