

Product datasheet

Normal Goat Serum (Sterile) ab138478

18 References

Overview

Product name	Normal Goat Serum (Sterile)
Host species	Goat
Tested applications	Suitable for: ICC/IF, IHC-P, IHC-Fr, Blocking
General notes	<p>Sterile goat serum for blocking in antibody-based applications or use in cell culture media. From mixed sex mixed strain goats. Origin: North America.</p> <p>Raised in: Goat</p> <p>Purity: Whole antiserum</p> <p>Sterilization: This product was aseptically filtered through a 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.</p>

Properties

Form	Liquid
Storage instructions	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Purity	Whole antiserum
Purification notes	<p>Sterilization: This product was aseptically filtered through a 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.</p>

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab138478 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.

Application	Abreviews	Notes
IHC-Fr		Use at an assay dependent concentration.
Blocking		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors