


## Product datasheet

# Ozagrel hydrochloride, thromboxane (TXA<sub>2</sub>) synthetase inhibitor ab142740

### Overview

<b>Product name</b>	Ozagrel hydrochloride, thromboxane (TXA <sub>2</sub> ) synthetase inhibitor
<b>Description</b>	Potent and selective thromboxane (TXA <sub>2</sub> ) synthetase inhibitor
<b>Biological description</b>	Potent and selective thromboxane (TXA <sub>2</sub> ) synthetase inhibitor (IC <sub>50</sub> = 4 nM). Shows no inhibitory activity at prostacyclin (PGI <sub>2</sub> ) synthetase, cyclooxygenase or PGE <sub>2</sub> isomerase (IC <sub>50</sub> > 1 mM). Inhibits platelet aggregation <i>in vitro</i> and arachidonate-induced arterial contraction <i>in vivo</i> .
<b>Purity</b>	> 98%
<b>CAS Number</b>	78712-43-3
<b>Chemical structure</b>	

### Properties

<b>Chemical name</b>	(E)-3-[4-(Imidazol-1-ylmethyl)phenyl]prop-2-enoic acid hydrochloride
<b>Molecular weight</b>	264.71
<b>Molecular formula</b>	C <sub>13</sub> H <sub>12</sub> N <sub>2</sub> O <sub>2</sub> .HCl
<b>PubChem identifier</b>	6438130
<b>Storage instructions</b>	Store at Room Temperature. The product can be stored for up to 12 months.
<b>Solubility overview</b>	Soluble in water to 100mM
<b>Handling</b>	<p>Wherever possible, you should prepare and use solutions on the same day. However, if you need to make up stock solutions in advance, we recommend that you store the solution as aliquots in tightly sealed vials at -20°C. Generally, these will be useable for up to one month. Before use, and prior to opening the vial we recommend that you allow your product to equilibrate to room temperature for at least 1 hour.</p> <p>Need more advice on solubility, usage and handling? Please visit our <a href="#">frequently asked questions (FAQ) page</a> for more details.</p>
<b>SMILES</b>	C1=CC(=CC=C1CN2C=CN=C2)C=CC(=O)O.Cl
<b>Source</b>	Synthetic

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES, NOT FOR USE IN HUMANS"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team