


Product datasheet

Anti-PADI3 / PAD3 antibody ab196288

1 Image

Overview

Product name	Anti-PADI3 / PAD3 antibody
Description	Rabbit polyclonal to PADI3 / PAD3
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat 
Immunogen	Recombinant full length protein within Human PADI3/ PAD3. The exact sequence is proprietary. Database link: BC109092
Positive control	Hela cell lysate and Human fetal skin lysate.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.4 Preservative: 0.02% Sodium azide Constituent: 99% PBS
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab196288** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000 - 1/2000. Predicted molecular weight: 75 kDa.

Target

Function	Catalyzes the deimination of arginine residues of proteins.
Tissue specificity	Hair follicles, and epidermis at very low levels.
Sequence similarities	Belongs to the protein arginine deiminase family.
Cellular localization	Cytoplasm.

Images



All lanes : Anti-PADI3 / PAD3 antibody (ab196288) at 1/1000 dilution

Lane 1 : HeLa cell lysate

Lane 2 : Human fetal skin lysate

Predicted band size: 75 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors